

DASEA® Ultramedia® Max mesenchymal stem cell serum-free medium



● INTRODUCTION

DASEA® Ultramedia® Max is the new generation product specifically developed for mesenchymal stem cell (MSC) culture. The product's composition is clearly defined and controllable, and the cultured cells exhibit good morphology and an ideal expansion rate, meeting the requirements for both primary and passaged culture needs.

* For research use only

● PRODUCT ADVANTAGE



Optimized Cell Morphology



Consistently High
Proliferation Rate



Desirable Cell Diameter

● PRODUCT SPECIFICATION

Model	Product Number	Description	Component Number	Package (per bottle)
Type A	RGM0071	Max basal medium (phenol red-free)	RGM1071	500mL
		Max supplement	RGM1073	15mL
Type B	RGM0072	Max basal medium (phenol red)	RGM1072	500mL
		Max supplement	RGM1073	15mL


● PERFORMANCE DATA

Passages	Seeding Density (cells/cm ²)	Period	Confluency	Harvested (cells) / T175	Expansion Fold	Total Harvested (cells)	Total Expansion Fold		
Primary	-	12-14 days	-	2.40E+06	-	1.68E+07	-		
P1	4000-5000	72h	80% ~ 90%	1.34E+07	15.37	2.58E+08	1.53E+01		
P2				1.24E+07	14.18	3.66E+09	2.17E+02		
P3				1.22E+07	13.96	5.11E+10	3.02E+03		
P4	5000-7000			1.41E+07	11.48	5.87E+11	3.47E+04		
P5				1.40E+07	11.46	6.72E+12	3.98E+05		
P6				1.29E+07	10.52	7.07E+13	4.18E+06		
P7				1.14E+07	9.3	6.58E+14	3.89E+07		
P8				8.92E+06	7.28	4.79E+15	2.83E+08		
P9				8.77E+06	7.16	3.43E+16	2.03E+09		
P10				8.99E+06	7.34	2.52E+17	1.49E+10		


The data mentioned refers to the average number of harvests and the proliferation multiples of six cell lines from umbilical cord, and the actual proliferation rate of the cell strains may be higher or lower than what is depicted in the figure.

● CELL PHENOTYPE

P3



P10



* Cell phenotype remains normal at P10



REGEN-αGEEK (Haining) Biotechnology Co., Ltd.

Tel: +860755-26412015

Web: <https://en.regengeek.com> Email: enquiry@regengeek.com

Address: No. 301, Building 6, Juanhu Science and Innovation Park, No. 500, Shuiyue Pavilion East Road, Xiashi Street, Haining City, Jiaxing City, Zhejiang Province, China.